

Regulations KD00887800

Disposal of Liquid Biosafety Waste

Disposal of liquid biosafety waste

The liquid waste will be autoclaved.

▲ Caution: do not add any kinds of chemical ingredients to the waste! No use of bleach! Bleach is corrosive and destroys the autoclave. In addition toxic vapors are generated

1. Before use label your waste bottles with your group name and floor. Mark your bottles with a max filling volume of **3 liters**. **Use the bottles to collect the liquid biosafety waste.**
2. Disconnect the waste bottle from your pump system when you reach the max. of 3 liters.
3. Close the bottle with the lid properly.
4. Pick up the trolley from the disposal room on your floor and bring the waste down to the V-Floor (LLF.F+6.014) in the waterproof box of the trolley.

Note: The decontamination cycle for liquids is limited to 3 liters of liquid. The decontamination cycle for larger volumes would take too long (more than 10 hours) and will not reach the necessary temperature for inactivation.

